

Widacherweg 7
Postfach 132
CH-8630 Rüti / ZH

Telefon 055 266 25 25
Tel. GSM 079 646 56 39
Telefax 055 266 25 26

E-Mail info@qelbrunner.ch
http:// www.qelbrunner.ch
http:// www.qelbrunner.com

Contained Powder Transfer Bag



As today's health and safety requirements are getting more stringent and demanding, the single-use powder transfer bag is the optimal solution to protect the operator from the Active Pharmaceutical Ingredients (API) and/or the API from operator-borne contamination.

The specific design of this flexible single-use transfer system allows for easy, safe and clean transfer of API powders and drug substances between different process steps.

The bottle shaped bag has a handle for easy transfer of powders and a PE tri-clamp connector for safe and easy connection.

Applications

- Transfer of API powders.
- Bulk Pharmaceutical Chemicals (BPC).
- Potent drug substances destined for clinical trials.
- Transfer of hazardous powders.
- Transfer to tablet or capsule machines.

Benefits

- Eliminates cleaning and validation steps.
- All contact materials are FDA approved.
- Liner materials comply with EC directive EEC 90/128, SI 1376, and FDA 21CFR.
- Low particle and bio-burden levels.
- 100% visual inspection.

Wir verweisen auf die technischen Datenblätter des Herstellers
Erstelldatum 11/03

Widacherweg 7
Postfach 132
CH-8630 Rüti / ZH

Telefon 055 266 25 25
Tel. GSM 079 646 56 39
Telefax 055 266 25 26

E-Mail info@qelbrunner.ch
http:// www.qelbrunner.ch
http:// www.qelbrunner.com

Contained Powder Transfer Bag



Features

- Cleanliness level 200 for interior and exterior of bags (according to IEST-STD-CC1246D).
- In-house PE film extrusion in a controlled environment.
- Bag is produced in an ISO Class 5 environment.
- Sterilizable by gamma-irradiation (25 kGray).
- Antistatic, pinhole-free LDPE film
- Tri-clamp connectors available in 2.5, 4 and 6 inches.
- Rigid handle for easy handling.

Wir verweisen auf die technischen Datenblätter des Herstellers
Erstelldatum 11/03