

Widacherweg 7  
Postfach 132  
CH-8630 Rüti / ZH

Telefon 055 266 25 25  
Tel. GSM 079 646 56 39  
Telefax 055 266 25 26

E-Mail qelbrunner@bluewin.ch  
http:// www.qelbrunner.ch  
http:// www.qelbrunner.com

## Klercide™ Isolator Cleaning Tool

A special autoclavable tool designed to take sterile disposable covers, which aids the cleaning of isolators and laminar flow cabinets.



### Advantages

#### Easy to use

Two part construction allows it to be easily passed through transfer hatches, and re-assembled. The tool is lightweight to aid reaching the rear of the cabinet and an extension handle is available for increased reach, especially for large or half suit isolators. Easier reaching helps alleviate back and upper limb disorders.

#### Efficient cleaning

The extra reach and specially shaped head ensure all areas of the Isolator, even difficult back corners can be cleaned reducing the risk of contaminants remaining. The design is suitable for cleaning flat surfaces, around fitted equipment and awkward corners.

#### Safe to use

The padded cover and rounded head ensures that it will not damage the internal surfaces of cabinets especially flexible film isolators. The Isolator Cleaning Tool is autoclavable and Klercide alcohols or biocides attack none of the component materials.

#### No residues

The materials used in the disposable pads do not contain any chemical binders so do not leave undesirable residues. The efficient cleaning helps reduce the risk of any disinfectant or spillage residues remaining.

Widacherweg 7  
Postfach 132  
CH-8630 Rüti / ZH

Telefon 055 266 25 25  
Tel. GSM 079 646 56 39  
Telefax 055 266 25 26

E-Mail qelbrunner@bluewin.ch  
http:// www.qelbrunner.ch  
http:// www.qelbrunner.com

**Description**

The Klercide Isolator Cleaning Tool comprises a two-part hard anodised aluminum handle and specially shaped head. They are cleaned with sterile alcohol and double bagged in cleanroom packaging before irradiation at no less than 25Kgy. The Isolator Cleaning Tool is fully autoclavable.

Once assembled the reach is approximately 600mm. A 300mm extension handle is available for larger cabinets or half suit isolators.

The Isolator Cleaning Tool is designed for use with sterile single use cleaning pads. The pads comprise of an inner core of absorbent viscose and a choice of a low particulate or 100% knitted polyester soft and strong cover. Each pad is supplied individually double bagged before irradiation at not less than 25Kgy.

Each new Klercide Isolator Cleaning Tool is supplied with 2 low particulate cleaning pads.

**Validation**

Manufactured by a BS EN ISO9001 accredited company.

Certificates of irradiation with every batch.

Shelf life of the cleaning pads is 3 years.

Further supporting information is available on this product please see the **Technical Information** section.

**Product Codes**

Code	Description	Unit of Sale
KCS3001	Klercide Isolator Cleaning Tool Supplied with 2 disposable cleaning pads KCS3055	1 unit
KCS3008	Klercide Isolator Cleaning Tool Extension Handle	1 unit
KCS3020	Klercide Sterile 100% Polyester Isolator Cleaning Pads Sterile, single use, 100% Polyester pads for use with KCS3001	20 pieces
KCS3055	Klercide Sterile Low Particulate Isolator Cleaning Pads Sterile, single use, Low Particulate pads for use with KCS3001	20 pieces

Code	Description	CHF Unit Price
KCS3001	Klercide Isolator Cleaning Tool Supplied with 2 disposable cleaning pads KCS3055	CHF 516.40 Pice
KCS3008	Klercide Isolator Cleaning Tool Extension Handle	CHF 166.80 Pice
KCS3020	Klercide Sterile 100% Polyester Isolator Cleaning Pads Sterile, single use, 100% Polyester pads for use with KCS3001	CHF 450.90 Case
KCS3055	Klercide Sterile Low Particulate Isolator Cleaning Pads Sterile, single use, Low Particulate pads for use with KCS3001	CHF 243.00 Case